产品说明书 Instructions for Use

CelliMax® MetaCHO™培养基 CelliMax® MetaCHO™ Medium

Cat. No.: CPDP050 (Powder)

Issued on August 1, 2024

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Introduction

- 适用于 CHO 细胞的培养
For CHO cell culture
- 化学成分确定
CD (Chemically Defined)
- 无动物源性成分
ADCF (Animal Derived Component Free)
- 不含血清
Serum free
- 不含酚红
Without Phenol Red
- 不含次黄嘌呤和胸苷
Without Hypoxanthine and Thymidine
- 不含 L-谷氨酰胺
Without L-Glutamine
- 不含水解物和重组蛋白
Without Hydrolysates and Recombinant Protein
- 含有 5.6g/L 葡萄糖
With 5.6g/L Glucose

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Storage

1. 干粉培养基密封保存在 2~8°C 的避光条件下。
Powder medium should be sealed at 2 to 8°C, away from light.
2. 配制成的液体培养基密封保存在 2~8°C 的避光条件下。
Liquid medium formulated from powder should be sealed dry at 2 to 8°C, away from light.

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Shelf Life*

1. 干粉培养基未开封情况下有效期 2 年, 开封后请尽快使用。详细信息请见标签。
Powder medium is valid for 2 years without opening. Please use as soon as possible after opening. Please refer to the label for details.
**上述有效期的制定, 均基于赛普生物的通用检测方法。客户可针对不同的使用用途或项目, 进行单独的有效期研究。
The shelf life is based on the general testing method of CelluPro. It is suggested to study on validity for different uses or projects.

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- ◆ 储液系统。
Liquid storage system.
- ◆ 除菌过滤系统。
Sterile filtration system.
- ◆ 为防止误差,不建议进行3升以下的小规模配制。
Small-scale preparation less than 3 liters is not recommended to prevent errors.
- ◆ 为确保产品质量,所有与培养基接触的器具或耗材,需保证无支原体和细菌内毒素污染。
All utensils or consumables in contact with the medium shall be mycoplasma and bacterial endotoxin contamination risk free for product quality.
- ◆ 与过滤后的液体培养基接触的过滤系统和储液系统需无菌。
Filtration system and liquid storage system in contact with the filtered liquid medium should be kept sterile.
- ◆ 建议增加 0.8 μ m 或 0.45 μ m 的预过滤,以加快过滤速度并提高过滤载量。
Pre filtration at 0.8 μ m or 0.45 μ m is recommended to improve filtration speed and filtration load.
- ◆ 搅拌溶解时间与配制规模和搅拌效率相关,本文中所有提及的搅拌时间仅为小规模配制时的参考值。若配制规模大或搅拌系统效率低,则需适当延长各搅拌步骤的时间,以确保可以达到充分的混合及溶解效果。
Dissolution time varies with preparation scales and mixing efficiency. All the stirring time herein is a reference for small-scale preparation. For larger preparation scale or lower mixing efficiency, each mixing time needs to extend appropriately to ensure sufficient mixing and dissolution effect.
- ◆ 过长的配制时间会带来更高的微生物负荷及细菌内毒素污染的风险,因此配制总时长不建议超过12小时。
Total preparation time is recommended not to exceed 12 hours, for long preparation time will pose higher risk of bioburden and endotoxin contamination.

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CelliMax® MetaCHO™ Medium is a protein-free CD basal medium. If the original medium with protein components (e.g. insulin and IGF) is replaced with CelliMax® MetaCHO™ Medium, corresponding substances can be added to shorten the cell adaptation time.

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*Except for process changes, other changes will not be notified.

©2024 CelluPro, Ltd. All rights reserved including data, graphics and images. CelliMax, CelluPro and Helix logo are registered trademarks of CelluPro, Ltd.